

Laser1 settings for FWHM scan

CURRENT (A)	Internal Attenuation (%)	Delay A (us)	Power Fast Acqiris	FWHM (ns) Slow	Center(ns) slow
20	97	7.2	4.95	39.7	11750
21	70	7.3	4.8	36.1	11750
21.5	58	7.335	4.9	34	11750
22	55	7.365	5	32.2	11750
22.5	52	7.387	5	31.2	11750
23	50	7.407	5	30.5	11750
23.5	46	7.418	5.2	29.1	11745
24	41	7.435	4.97	27.7	11750
24.5	34	7.45	4.93	26.6	11750
25	30	7.46	5.2	25.6	11750
25.5	28	7.47	5.1	24.7	11750

Reference done at 23A, start working point.

Measurement of power, pulse FWHM and Center done after very short run.

22 Sept 2009